



ICVR 2024

2024 the 10th International Conference on Virtual Reality



Bournemouth, United Kingdom
July 24-26, 2024

A VISUO-HAPTIC FEEDBACK SURGICAL SIMULATOR FOR TWIN-TO-TWIN TRANSFUSION SYNDROME




**MICHAEL KASMAN, TRISTAN ALKIS, DR. KENNETH J. MOISE JR., DR. MICHAEL BEBBINGTON,
AND PROF. ANN MAJEWICZ FEY**

TEXAS Robotics




The University of Texas at Austin
Dell Medical School



What is Twin-to-Twin Transfusion Syndrome (TTTS)?

- Approximately 10-15% of monochorionic twin pregnancies
- Life-threatening pregnancy complication where there is unequal distribution of blood between the twins
- The recipient twin develops excessively, and the donor is at high risk of cardiac failure
- Left untreated, can lead to the death of one or both twins
- Rare condition, incredibly difficult procedure, and hard to train



Fetoscopic Laser Photocoagulation (FLP)

The **goal** is to **identify anastomoses** and **ablate vessels** to **decouple the twins' vascular system**

1. Final all anastomosis on the surface of the placenta
2. Ablate each anastomosis site
3. “Solomon method” – create a laser path connecting each ablation site

1



2




3




Improving Proprioception

Where are the surgeons in the workspace relative to the identified anastomoses?

Real Fetoscopic Procedure




Placenta Anatomical Features



- Challenging to maintain an accurate sense of position in relation to identified anastomoses
- Difficult to remember vein and artery patterns

Type of TTTS Simulators

Virtual



ADVANTAGES

- Simulates the baby and placenta
- Custom fetoscope interface attached to an abdomen model
- Fetoscope settings (brightness, laser power, etc.)
- UI/UX

DISADVANTAGES


- Absence of camera noise, unrealistic textures, laser effect lacks differentiation between vessels and tissue
- Placental vessels are flat
- No haptic or visual feedback

Type of TTTS Simulators

ADVANTAGES

- Feels realistic and functional
- Can render placental image for anterior and posterior locations
- Able to add visual indicators to the placenta
- EM tracker at tip for motion analysis

Mixed-Reality



DISADVANTAGES

- Does not render essential elements: twins, fluid, particles, laser effects, etc.
- Placenta is a planar rendering
- No haptic feedback

Type of TTTS Simulators


ADVANTAGES

- High fidelity
- Placental vasculature injected with dye to differentiate between veins and arteries

DISADVANTAGES


- Utility is post-delivery, limited placenta availability
- Does not simulate floating twins and amniotic particles
- Very expensive

Real Environment



Contributions of This Work

1. The development of a hyper-realistic simulator with all essential elements of TTTS that fetal surgeons will use
2. The improvement of surgical education with visuo-haptic feedback and assistance
3. Proof of concept of assistance features for real surgery



No Assistance

Visual Guidance

Haptic Guidance

User

Target



$$k * \min_{1 \leq j \leq m} |\vec{x} - \vec{x}_j|$$

Mini-Map Guidance

Combined Guidance

HoloLens Guidance

HoloLens Trainer Demo

A person wearing a HoloLens headset is standing at a desk in a training environment. They are interacting with a laptop and a monitor. A hand is visible in the foreground, holding a small device. The desk has various items, including a mug, a bowl, and a small sign. The background shows a desk with a monitor and a chair.

Experimental Setup


- **Participants:** 30 novices, 1 expert surgeon

- 5 novices for each assistance case

- **Age (years old)**

- **Expert:** 68
 - **No Assistance:** 23 ± 5
 - **Haptic Guidance:** 32 ± 17
 - **Visual Guidance:** 32 ± 17
 - **Mini-Map Guidance:** 26 ± 4
 - **Combined Guidance:** 24 ± 4
 - **HoloLens Guidance:** 24 ± 4


- Anastomosis waypoints and trajectory of laser ablation task are collected for analysis




Experimental Task

Participants repeated these tasks 2 times: **first trial (no assistance)** and **second trial (assistance)**


1. Identify anastomosis sites
2. Laser ablation of anastomosis sites
3. Create a laser path



(a) Optimal Placenta Path



(b) Reversed Placenta Path



Experimental Design

Pre-Experiment Survey



Surgical Video of TTTS Procedure



Instruction of Task



Testing/Training Phase



First and Second Trial



Post-Experiment Survey

Metrics

1. Waypoints Placed

$$\sum_{i=1}^k \sum_{j=1}^n |\vec{x}_i - \vec{x}_j|^2 \quad (\text{Eq. 1})$$

$$(x_i - x_{i+1})/x_i < 0.2 \quad (\text{Eq. 2})$$

2. Laser Path Deviation

$$\frac{1}{n} \sum_{i=1}^n \min_{1 \leq j \leq m} |\vec{x}_i - \vec{x}_j| \quad (\text{Eq. 3})$$

3. Velocity, Acceleration, Jerk

$$v_{\text{avg}} = \frac{1}{n} \sum_{i=1}^{n-1} \frac{|\vec{x}_{i+1} - \vec{x}_i|}{\Delta t_i} \quad (\text{Eq. 4})$$

4. Efficiency

$$\frac{1}{\text{Task Completion Time}} \quad (\text{Eq. 5})$$


5. Waypoint and Laser Error

Same computation as Eq. 3 except comparing using **the novices' waypoints to the expert's waypoints** and **the novices' laser ablation path to the expert's path**.


Results

No Assistance


First Trial




Second Trial





First Trial





Second Trial





Haptics





Mini-Map



Combined Assistance



HoloLens




Result – Difference Between Groups

- Kruskal-Wallis test** does not show a significant difference between groups
- However, individual p-value compared between specific groups show statistical significance using **post-hoc Dunn test**
 - Number of Waypoints Placed
 - Laser Path Deviation
 - Waypoint Error
 - Laser Error

Metrics	Session		Guidance Cases	
	p	post-hoc comparisons	p	post-hoc comparisons
Waypoints Placed	0.926	(Trial One < Trial Two)	0.133	(C < V < MM < NA) < H < (HL)
Laser Path Deviation	0.345	(Trial Two < Trial One)	0.202	HL < (MM < C < V) < NA < H
Waypoint Error	0.449	(Trial Two < Trial One)	0.100	HL < C < (MM < H = NA) < V
Laser Error	0.092	(Trial Two < Trial One)	0.184	HL < (V < C) < NA < (H) < MM


Asterisk * represents significant differences between groups $p < 0.05$. Round bracket denotes no significance between the pairs of group levels.



Result – Task Efficiency and Velocity

- Efficiency and Velocity: Significant differences between trial one and trial two; task is becoming easier to do

Metrics	Session		Guidance Cases	
	<i>p</i>	post-hoc comparisons	<i>p</i>	post-hoc comparisons
Velocity	0.021*	Trial One < Trial Two	0.758	(V < MM < C < HL < NA < H)
Efficiency	0.038*	Trial One < Trial Two	0.701	(V < C < MM < NA < H < HL)



User Feedback

NASA
TLX

		E	N	V	MM	HL	H	C
Q1:	How mentally demanding was the task?	65	60 (30)	86 (21)	74 (22)	58 (10)	70 (6)	70 (20)
Q2:	How physically demanding was the task?	20	59 (39)	21 (29)	60 (26)	30 (18)	33 (19)	31 (22)
Q3:	How hurried or rushed was the pace of the task?	10	28 (26)	23 (25)	35 (16)	33 (19)	37 (18)	17 (8)
Q4:	How successful were you in accomplishing what you were asked to do?	60	62 (20)	56 (33)	44 (18)	37 (13)	65 (13)	43 (22)
Q5:	How hard did you have to work to accomplish your level of performance?	45	64 (21)	83 (21)	72 (13)	52 (15)	64 (19)	72 (10)
Q6:	How insecure, discouraged, irritated, stressed, and annoyed were you?	25	51 (35)	60 (31)	51 (35)	33 (24)	58 (21)	41 (39)
Q7:	The simulator looks realistic	80	76 (17)	80 (16)	85 (10)	75 (10)	76 (17)	84 (17)
Q8:	The simulator feels realistic	60	80 (20)	88 (11)	84 (9)	85 (10)	84 (17)	76 (22)
Q9:	The content of the simulation is relevant to the surgical procedure	60	92 (11)	100 (0)	85 (19)	85 (19)	80 (14)	96 (9)
Q10:	The simulation was fun to complete	100	76 (33)	84 (17)	64 (17)	95 (10)	76 (22)	68 (30)
Q11:	I was anxious during the task	100	56 (41)	60 (20)	48 (23)	55 (34)	52 (27)	84 (17)
Q12:	The task was difficult to complete	30	80 (24)	76 (22)	84 (9)	75 (10)	80 (14)	68 (18)
Q13:	Age	68	23 (5)	32 (17)	26 (4)	24 (4)	32 (17)	24 (4)
Q14:	Experience with human/machine interactive devices	1	2.8 (1.1)	1.8 (1.3)	2.4 (0.9)	2.8 (1.1)	2.2 (1.1)	3.2 (0.9)

Data is shown as mean (standard deviation). Q1-6 are NASA TLX questions. For Q1-4,6: 0 is very low and 100 is very high. For Q5: 0 is perfect and 100 is failure. For Q7-10,12: 0 is strongly disagree and 100 is strongly agree. For Q11: 100 is not anxious at all while 0 is very anxious. For Q14: 0 is none and 4 extensive. E=Surgeon. Assistance cases: N=None, V=Visual, MM=Minimap, HL=HoloLens, H=Haptic, C=Combined.

Pre-Exp
Survey

Expert Surgeon Feedback

TABLE II: Post-Experiment Survey: Surgeon-Specific Questions

Q1:	The simulator resembles real-life surgery	4
Q2:	The simulator can differentiate between different performance levels (e.g. novice vs expert)	4
Q3:	The simulator can accurately assess the skill level of the user	4
Q4:	This simulator is better than others or expected	4
Q5:	This simulator gives results consistent with other tests for the same skill	4
Q6:	The performance in the simulator would translate to a similar performance in a real-life procedure	3
Q7:	A surgeon who performs well in the simulation will perform well in a real surgery	3
Q8:	Another evaluator would assess a surgeon's performance on the simulator similar to yours	4
Q9:	If you were to repeat the test at a different time you would achieve a similar result	5
Q10:	This simulator is good for learning	4
Q11:	I would recommend this simulator to a peer	5
Q12:	Confirm position of placenta, fetuses, and cord insertions	3
Q13:	Identification of intertwin dividing membrane	3
Q14:	Mapping of placental surface and vascular equator	4
Q15:	Use of selective laser techniques (only coagulation of anastomoses)	4
Q16:	Use of Solomon technique	4
Q17:	Identify and record number and type of anastomosis sites coagulated	5

1-11 are rated strongly agree or strongly disagree, 5 and 0 respectively. 12-17 are rated as how accurate, from perfect to not at all, 5 and 0 respectively.

Limitations and Future Work

- Lack of situational haptics for user interface
- Overall improve laser depth perception
- Combining a physical model with the virtual simulation
- Mimic the energy loss when the laser is not angled at approximately 0°
- Intrauterine bleeding onto the placenta
- A viscous field between points that directs you back to the path that connects the identified anastomoses

Conclusion

- In this work we introduce a novel hyper-realistic TTTS simulator
- Through a human subject study, our simulator could improve performance of laser ablation path accuracy
- Significant difference between certain assistance cases cases
- Significant difference between task efficiency and laser ablation velocity from trial one to trial two
- The results of our simulator could improve surgical training and inform the development of novel robotic systems for TTTS

Thank you!
Any Questions?